



IHE Integration Statement

PaxeraServer

2020



INTRODUCTION

Overview

The Integrating the Healthcare Enterprise (IHE) Integration Statement was developed as a high level statement of integration for products, which were intended to interoperate, based on IHE Integration Profiles.

The IHE Technical Framework defines a subset of the functional components of the healthcare enterprise, called IHE Actors, and specifies their interactions in terms of a set of coordinated, standards-based transactions. These transactions are organized into functional units called Integration Profiles that highlight their capacity to address specific clinical needs.

This document identifies the IHE Actors, IHE Integration Profiles and options PaxeraHealth has implemented in each product version.

Users familiar with these concepts can use Integration Statements as an aid to determine what level of integration this product supports with other systems and what clinical and operational benefits such integration might provide.

For further exploration, additional information can be found in the DICOM and HL7 Conformance Statements of this product and in the IHE Technical Framework.

Important Note to the Reader

This IHE Integration Statement does not guarantee successful interoperability of this product with other vendor's products. It is the user's responsibility to thoroughly analyze the application requirements and to specify a solution that integrates the PaxeraServer with other vendor's software and systems.

The information presented in this document is correct at the time of going to press. However, it is subject to change without prior notice. This document may not be photocopied, reproduced or translated into another language, in whole or in part, without the consent of Paxerahealth Corp.

Future Evolution

Paxerahealth understands that the IHE Technical Framework and the underlying communication standards will evolve to meet the user's growing requirements.

Paxerahealth is actively involved in this evolution and plans to adapt its products to future versions of the IHE Technical Framework. Such evolutions may require changes to already installed systems



IHE INTEGRATION STATEMENT

PaxeraHealth IHE Integration Statement				
Vendor	Product Name	Version	Date	
PaxeraHealth Corp.	PaxeraServer PaxeraUltima PaxeraVNA	7.7.x	Feb 05 2020	

This product implements all transactions required in the IHE Technical Framework to support the IHE profiles, actors and options listed below.

Profiles Implemented	Actors Implemented	Options Implemented	
	RAD – Radiology Domain		
Digital Breast Tomosynthesis (DBT)	Image Manager/Image Archive	None	
	Image Display	None	
Imaging Object Change Management (IOCM)	Change Requester	None	
	Image Manager/Image Archive	None	
	Image Display	None	
Portable Data for Imaging (PDI)	Portable Media Creator	DVD Media	
Cross-enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer	None	
	Imaging Document Source	Set of DICOM Instances	
Mammography Image (MAMMO)	Image Manager/Image Archive	None	
	Image Display	None	
Scheduled Workflow (SWF)	Image Manager/Image Archive	None	
	Performed Procedure Step Manager	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive	None	
	Performed Procedure Step Manager	None	



Key Image Note (KIN)	Image Manager/Image Archive	None	
	Image Display	None	
	Evidence Creator	None	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive	None	
	Image Display	None	
	Evidence Creator	None	
Access to Radiology Information(ARI)	Image Manager/Image Archive	None	
ІТ	I – IT Infrastructure Domain		
Audit Trail and Node Authentication (ATNA)*	Secure Application	Radiology Audit Trail	
Consistent Time (CT)	Time Client	None	
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	None	
	dress for vendor's IHE information ocs/IHE Integration Statement for		
Links to Standards Co	nformance Statements for the	Implementation	
	https://www.paxerahealth.com/docs/HL7_doc.pdf		
HL7			

 $[\]ensuremath{^{\boldsymbol{*}}}$ Limitted availability, Only available for customer upon request for this version.



PaxeraHaelth Corp. Contacts

300 Washington St, Newton, MA 02458, USA

Tel: 1-888-430-0005 Fax: 1-816-817-4648

www.paxerahealth.com

Sales: sales@paxerahealth.com
support@paxerahealth.com